

PAJUNK®

*NRFit™ - New ISO standard
Products and Components
for Neuraxial Anesthesia*



NRFit™

Enhancing Patient Safety

Made in
GERMANY

NRFit™ – New ISO Standard

Enhanced patient safety in Neuraxial Anesthesia

Due to various incidents of misconnections on various supply lines to the patient, the new ISO standard 80369-6 has been developed and approved internationally in 2016. This standard defines new requirements for small bore connectors in the field of Regional Anesthesia and intends to prevent the risk of accidental misconnections and injections, thus improving patient safety. The new connector following the ISO 80369-6 standard is termed NRFit™.

Commencing January 1st, 2017, California law prohibits a health facility from using an Epidural connector that would fit into a connector other than the type it was intended for, unless an emergency or urgent situation exists and the prohibition would impair the ability to provide health care. **This means California health facilities should no longer use epidural products with Luer connectors as of January 1st, 2017.**

As a leading manufacturer in the field of Regional Anesthesia, PAJUNK® sees the necessity to advance these measures to improve patient safety and therefore will provide a full portfolio of NRFit™ products for Spinal, Epidural and CSE procedures in 2017. Products for Peripheral Nerve Blocks and Continuous Wound Infusion will be available at a later date.


PAJUNK® will actively accompany and support the transition process in your facility! Our staff is ready to advise and assist you with upcoming NRFit™ product questions!



SPROTTE® - Spinal fit with NRFit™ connectors



EpiLong Soft - Epidural fit with NRFit™ connectors



New HUB connectors:

All Regional Anesthesia products will be equipped with the new NRFit™ connectors.

➡ The procedure remains the same, but patient safety is being enhanced.



Color Coding:

The components color code for the NRFit™ range of products will be yellow.

➡ This allows for easy identification of the correct connector types.



Package Labeling:

All individual and secondary packaging of NRFit™ products will come with a clear, easy to read labeling.

➡ Additional labeling helps to prevent mix-ups - especially in the early stage of the conversion process.

NRFit is a registered trademark of GEDSA used with their permission.

PAJUNK® MEDICAL SYSTEMS L.P. • 6611 Bay Circle, Suite 140 • Norcross, GA 30071, USA
Phone: +1.888.9PAJUNK (72.5865) • Fax +1.678-514-3388 • E-Mail: info@pajunk-usa.com