


PAJUNK®

Paediatric Anaesthesia

*Cannula systems for caudal,
epidural, and plexus anaesthesia*



Paediatric anaesthesia



Fine and precise

The PAJUNK® Special Program for Paediatric

As manufactured of various specialty cannula, PAJUNK® offers a complete product line designed specifically for paediatric anaesthesia in

regional anaesthesia. We offer high quality cannula for caudal, spinal, epidural, and plexus anaesthesia.

Caudal Anaesthesia

SPROTTE® 27 G



The SPROTTE® has initiated a worldwide trend to atraumatic puncture, which has proven itself in paediatric anaesthesia applications as well.

Crawford cannula 23 G, graduated, with depth stop




The Caudal cannula with its cone point geometry according to Crawford has proven to be particularly suitable in paediatric puncture.




Epidural Anaesthesia

PAJUNK® produces two special types of cannulas for Epidural anaesthesia:


With the well proven Tuohy tip:



EpiLong Paed with the well proven Tuohy tip:




SPROTTE® cannula



Crawford tip



Tuohy tip




Tuohy tip


ric Anaesthesia

Plexus Anaesthesia


UniPlex NanoLine cannula „Single Shot“: (UP/UPS) cannula




The UniPlex NanoLine cannulas were developed by PAJUNK® specifically for the disposable application in the peripheral conductive anaesthesia. These cannulas come in many different gages and sizes. PAJUNK® understands the growing trend of this application in paediatric anaesthesia so we offer two versions a SPROTTE® tip and a Facet tip with a shorter lengths and very small gages.



SPROTTE® tip




Facet tip



SPROTTE® SPECIAL tip

PlexoLong NanoLine according to Meier „Catheter Technique“



The patented Plexus Catheter Set was developed in cooperation with Dr. Meier specifically for the continuous peripheral techniques for regional anaesthesia for children.



Caudal

Product	Size	Item No.	PU
SPROTTE® cannula	27 G x 35 mm	111151-27A	25
SPROTTE® cannula	25 G x 35 mm	001151-29E	25
SPROTTE® cannula	24 G x 35 mm	001151-30E	25
SPROTTE® cannula	22 G x 50 mm	071151-30C	25
Crawford cannula	23 G x 25 mm	1149-12D025	25
Crawford cannula	22 G x 40 mm	1149-12E040	25

Epidural

Product	Size	Catheter	X-ray contrast catheter Guide 0,2 µm flat filter Tuohy Borst adapter LOR syringe	Item No.	PU
EpiLong Set I					
Tuohy cannula	19 G x 45 mm	23 G x 50 cm	• • • •	0131153-50	10
Epidural cannula					
Tuohy cannula	19 G x 50 mm			1150-4I050	25
Tuohy cannula	18 G x 50 mm			1150-4K050	25
EpiLong Paed					
Tuohy NanoLine cannula	21 G x 40 mm	25 G x 50 cm	• • • •	0131156-51	10

Plexus

Product	Size	Catheter	Item No.	PU
PlexoLong Set				
with SPROTTE® SPECIAL tip	21 G x 30 mm	25 G x 50 cm	541156-31B	10
UniPlex NanoLine Facet Tip				
UP 3/25	24 G x 25 mm		001156-75	10
UP 3/40	22 G x 40 mm		001156-70	
UP 3/50	22 G x 50 mm		001156-74	10
UniPlex NanoLine SPROTTE®				
UPS 3/40	24 G x 40 mm		001156-30G	10
UPS 3/50	25 G x 50 mm		001156-29G	10

PAJUNK Medical Systems L.P.
 6611 Bay Circle, Suite 140
 Norcross, Georgia 30071
 Phone: 770.493.6832
 Fax: 678.514.3388
 E-mail: info@pajunk-usa.com