

Regional Anesthesia Procedure Trays

PAJUNK trays provide quality and convenience for your procedures.



Details Separate the Good from the Great

As a market leader in Regional Anesthesia, PAJUNK® offers a full range of procedural tray solutions for Spinals, Epidurals, CSEs and Nerve Blocks. Building on feedback from clinical users, we've created a standard portfolio of the most popular tray configurations on the market. If you need something different, you have the CHOICE to create your own Custom Tray to fit your procedural needs.



Small Footprint Trays to Save Storage Space


All Pajunk full trays come in a compact design to allow for the smallest possible footprint. Comparable competitive full trays use up to 30% more shelving space.

- PAJUNK TRAYS: 10" x 10"
- COMPETITOR 1: 10.5" x 12" - 26% Bigger
- COMPETITOR 2: 10" x 13" - 30% Bigger



4" (103mm) Procedure Needles

Pajunk provides longer-length Neuraxial needles as a standard to eliminate the need for dropping on longer needles, and to meet the procedural needs of most patients, regardless of size. This also improves patient comfort by avoiding hubbing the needle when possible.



Easy to Read Tray Labels for Fast and Easy Identification

Color-coded labels simplify the identification process when pulling trays from shelves. All major tray information is easy to read and componentry is grouped for ease of reference to quickly identify which procedure tray it is.

- Preparation
- Accessories
- Procedure
- Medication



Original Pajunk Components

High-Quality Pajunk Components with Unique Features

- **Spinal** SPROTTE®
- **Epidural** EpiLong/EpiLong Soft
- **CSE** EpiSpin Lock
- **PNB/Single Shot** SonoPlex II
- **CPNB** SonoLong Catheters (CTN)
- **CPNB** E-Cath/E-Cath Plus (CON)

PAJUNK® CHOICE Procedure Trays

PAJUNK® offers a Tray Inventory Guarantee Program (TIP) to ensure you get the products you need for your procedures **WHEN** you need them. Ask your sales representative how you can get started with guaranteed inventory!



Support Trays

Pajunk offers a wide range of products to support your procedure if you're using various sizes of procedural needles.



Custom Trays

Let us know what Pajunk products and supporting components you need to create the perfect tray solution for you.

Spinal Trays

Together, with the founders of Pajunk, Prof. Dr. Sprotte developed the SPROTTE® atraumatic spinal needle, which revolutionized and revived the field of Spinal Anesthesia in the early 1980's, allowing clinicians to perform safer spinals by significantly reducing PDPH rates.

→ Sprotte

Epidural Trays

Pajunk offers a wide range of needles and catheters for safe Epidural Anesthesia, including Tuohy needles in different sizes and lengths, as well as EpiLong Nylon and EpiLong Soft Springwound catheters.

Pajunk Tuohy needles are made of high-quality polished steel, making them the perfect fit for Epidural Anesthesia.

- EpiLong Soft
- Tuohy

Nerve Block and CPNB Trays

Whether a standard tray, or one you create yourself, you'll be getting the best single shot and continuous nerve block products on the market. Pajunk is well known for highly echogenic, dual-guided innovations under ultrasound for both single-shot and continuous nerve block products: SonoPlex II and E-Cath/E-Cath Plus

- SonoPlex® II
- SonoLong Catheters (CTN)
- E-Cath®/E-Cath Plus (CON)



Combined Spinal Epidural (CSE) Trays

With EpiSpin Lock, Pajunk combines two techniques into one easy-to-use product: Spinal Anesthesia using the atraumatic Sprotte needle, and continuous epidural catheter placement using the Tuohy needle. The EpiSpin Lock offers smooth spinal puncture while minimizing needle tip damage through the back-eye opening on the distal end of the Tuohy needle. With the integrated locking mechanism, the Sprotte needle is stabilized for the injection of anesthetic while reducing the risk of needle displacement.

- EpiSpin Lock
- Sprotte

Procedure Trays

All trays come with the displayed components; the Pajunk needle included in each tray is listed under the specific item numbers in the Table. *All tray configurations, including custom, are available in NRFit®

Single Shot Nerve Block Trays

Components

3ml Tinted Chloraprep™ Applicator, 18Gx1.5" Hypodermic Needle, 25Gx1.5" Hypodermic Needle, 3 Gauze Pads, Paper Towel, CSR Wrap, 5ml Luer Lock Syringe, 20ml Luer Lock Syringe, Needle Foam Block, Drug ID Label, 5ml 1% Lidocaine HCl

Luer Item No.	NRFit Item No.	Tray Description
TDL101	TDN101	SonoPlex II 21G x 100mm (sterile item code 001287-77)
TDL131	TDN131	SonoPlex II 21G x 100mm (sterile item code 001287-77) with 11" (28cm) Safersonic probe cover and ultrasound gel
TDL102	TDN102	SonoPlex II 22G x 50mm (sterile item code 001287-74)
TDL132	TDN132	SonoPlex II 22G x 50mm (sterile item code 001287-74) with 11" (28cm) Safersonic probe cover and ultrasound gel

Single Shot Nerve Block Support Trays

Components

3ml Tinted Chloraprep™ Applicator, 18Gx1.5" Hypodermic Needle, 25Gx1.5" Hypodermic Needle, 3 Gauze Pads, Paper Towel, CSR Wrap, 5ml Luer Lock Syringe, 20ml Luer Lock Syringe, Needle Foam Block, Drug ID Label, 5ml 1% Lidocaine HCl

Luer Item No.	NRFit Item No.	Tray Description
TDL109	TDN109	Single Shot Support Tray, No Pajunk Components
TDL139	TDN139	Single Shot Support Tray, No Pajunk Components with 11" (28cm) Safersonic probe cover and ultrasound gel

Continuous Nerve Block Support Trays

Components

3ml Tinted Chloraprep™ Applicator, 18Gx1.5" Hypodermic Needle, 25Gx1.5" Hypodermic Needle, 3 Gauze Pads, Paper Towel, CSR Wrap, 4" Filterstraw, 3ml Luer Lock Syringe, 5ml Luer Lock Syringe, 20ml Luer Lock Syringe, 24"x32" Clear, Perforated Drape with 4"x6" Rectangular Opening, 2 4"x4.5" Film Dressings, 6 Wound Closure Strips, Skin Protectant, Needle Foam Block, Drug ID Label, 5ml 1% Lidocaine HCl

Luer Item No.	NRFit Item No.	Tray Description
TDL110	TDN110	Continuous Support Tray, No Pajunk Components
TDL140	TDN140	Continuous Support Tray, No Pajunk Components with 47" (120cm) Safersonic probe cover and ultrasound gel



Continuous Nerve Block Trays with E-Cath/Plus Catheter-Over-Needle (CON)

Components

3ml Tinted Chloraprep™ Applicator, 18Gx1.5" Hypodermic Needle, 25Gx1.5" Hypodermic Needle, 3 Gauze Pads, Paper Towel, CSR Wrap, 4" Filterstraw, 3ml Luer Lock Syringe, 5ml Luer Lock Syringe, 20ml Luer Lock Syringe, 24"x32" Clear Drape with 4"x6" Rectangular Opening, 2 4"x4.5" Film Dressings, 3 Wound Closure Strips, FixoCath Catheter Fixation Device, Skin Protectant, Needle Foam Block, Drug ID Label, 5ml 1% Lidocaine HCl

Luer Item No.	NRFit Item No.	Tray Description
TDL103	TDN103	E-Cath 18G x 51mm (2") indwelling catheter, SonoPlex 21G x 68mm (2-2/3") (sterile item code 201285-40E)
TDL133	TDN133	E-Cath 18G x 51mm (2") indwelling catheter, SonoPlex 21G x 68mm (2-2/3") (sterile item code 201285-40E) with 47" (120cm) Safersonic probe cover and ultrasound gel
TDL104	TDN104	E-Cath 18G x 83mm (3-1/4") indwelling catheter, SonoPlex 21G x 100mm (4") (sterile item code 241285-40E)
TDL134	TDN134	E-Cath 18G x 83mm (3-1/4") indwelling catheter, SonoPlex 21G x 100mm (4") (sterile item code 241285-40E) with 47" (120cm) Safersonic probe cover and ultrasound gel
TDL111	TDN111	E-Cath Plus 18G x 51mm (2") indwelling catheter, SonoPlex 21G x 68mm (2-2/3") (sterile item code 201285-41E)
TDL141	TDN141	E-Cath Plus 18G x 51mm (2") indwelling catheter, SonoPlex 21G x 68mm (2-2/3") (sterile item code 201285-41E) with 47" (120cm) Safersonic probe cover and ultrasound gel
TDL112	TDN112	E-Cath Plus 18G x 83mm (3-1/4") indwelling catheter, SonoPlex 21G x 100mm (4") (sterile item code 241285-41E)
TDL142	TDN142	E-Cath Plus 18G x 83mm (3-1/4") indwelling catheter, SonoPlex 21G x 100mm (4") (sterile item code 241285-41E) with 47" (120cm) Safersonic probe cover and ultrasound gel
TDL135	TDN135	E-Cath 18G x 100mm (4") indwelling catheter, SonoPlex 21G x 118mm (4-2/3") (sterile item code 261285-40E) with 47" (120cm) Safersonic probe cover and ultrasound gel
TDL145	TDN145	E-Cath Plus 18G x 100mm (4") indwelling catheter, SonoPlex 21G x 118mm (4-2/3") (sterile item code 261285-41E) with 47" (120cm) Safersonic probe cover and ultrasound gel (Clear Drape is Perforated)

Continuous Nerve Block Trays with SonoLong Catheter-Through-Needle (CTN)

Components

3ml Tinted Chloraprep™ Applicator, 18Gx1.5" Hypodermic Needle, 25Gx1.5" Hypodermic Needle, 3 Gauze Pads, Paper Towel, CSR Wrap, 4" Filterstraw, 3ml Luer Lock Syringe, 5ml Luer Lock Syringe, 20ml Luer Lock Syringe, 24"x32" Clear Drape with 4"x6" Rectangular Opening, 2 4"x4.5" Film Dressings, 6 Wound Closure Strips, Skin Protectant, Needle Foam Block, Drug ID Label, 5ml 1% Lidocaine HCl

Luer Item No.	NRFit Item No.	Tray Description
TDL105	TDN105	SonoLong Echo Tuohy 18G x 100mm 60cm injection tubing, 60cm stimulation cable, SonoLong Echo 20G x 23" (60cm) echogenic catheter with closed tip and 3 lateral openings, catheter container & introducer aid, and filter with clamping adapter (sterile item code 521185-31CU)
TDL106	TDN106	SonoLong Echo Tuohy 18G x 50mm 60cm injection tubing, 60cm stimulation cable, SonoLong Echo 20G x 23" (60cm) echogenic catheter with closed tip and 3 lateral openings, catheter container & introducer aid, and filter with clamping adapter (sterile item code 531185-31CU)
TDL107	TDN107	SonoLong Sono Tuohy 18G x 100mm 60cm injection tubing, 60cm stimulation cable, SonoLong Sono 20G x 20" (50cm) echogenic catheter with closed tip, helical coil, steel stylet, catheter container & introducer aid, and filter with clamping adapter

Procedure Trays

All trays come with the displayed components; the Pajunk product included in each tray is listed under the specific item numbers in the Table. *All tray configurations, including custom, are available in NRFit®

Spinal Trays

Components

25Gx1.5" Hypodermic Needle, 18Gx1.5" Hypodermic Needle, 1.75" Filterstraw, 3ml Luer Lock Syringe, 5ml Luer Slip Syringe, 3/4 fl. oz. Povidone-Iodine, 3 Prep Sponges, 24"x32" Clear Drape with 4" Circular Opening, 3 Gauze Pads, Paper Towel, CSR Wrap, Needle Foam Block, Drug ID Label, 5ml 1% Lidocaine HCl, 2ml 0.75% Bupivacaine HCL with 8.25% Dextrose OR Marcaine, where indicated

Luer Item No.	NRFit Item No.	Tray Description
TBL121	TBN121	24G x 103mm SPROTTE, 40mm introducer, DRY (No Spinal Drug) (sterile item code 021151-30A)
TBL131	TBN131	24G x 103mm SPROTTE, 40mm introducer, with Spinal Drug (sterile item code 021151-30A)
TBL151	TBN151	24G x 90mm SPROTTE, 30mm introducer, with Spinal Drug (sterile item code 121151-30A)
TBL122	TBN122	25G x 103mm SPROTTE, 40mm introducer, DRY (No Spinal Drug) (sterile item code 041151-29A)
TBL132	TBN132	25G x 103mm SPROTTE, 40mm introducer, with Spinal Drug (sterile item code 041151-29A)
TBL152	TBN152	25G x 90mm SPROTTE, 30mm introducer, with Spinal Drug (sterile item code 021151-29A)
TBL134	TBN134	27G x 103mm SPROTTE, 40mm introducer, with Spinal Drug (sterile item code 141151-27A)
TBL154	TBN154	27G x 90mm SPROTTE, 30mm introducer, with Spinal Drug (sterile item code 121151-27A)



Continuous Epidural Trays

Components

Clamping Adapter, Epidural Filter, Epidural Label, 25Gx1.5" Hypodermic Needle, 18Gx1.5" Hypodermic Needle, 4" Filterstraw, Glass LOR 5ml Luer Lock Syringe, 3ml Luer Lock Syringe, 20ml Luer Lock Syringe, 3/4 fl. oz. Povidone-Iodine, 3 Prep Sponges, 4"x4.5" Film Dressing, 24"x32" Clear Drape with 4" Circular Opening, 3 Gauze Pads, Paper Towel, CSR Wrap, Needle Foam Block, Drug ID Label, 5ml 1% Lidocaine HCl, 5ml Lidocaine+Epinephrine (Test Dose), 10ml Sodium Chloride

Luer Item No.	NRFit Item No.	Tray Description
TAL901	TAN901	17G x 103mm Tuohy, 19G Soft Closed Tip Catheter, DRY (No Test Dose)
TAL101	TAN101	17G x 103mm Tuohy, 19G Soft Closed Tip Catheter
TAL103	TAN103	18G x 103mm Tuohy, 20G Nylon Closed Tip Catheter
TAL151	TAN151	17G x 90mm Tuohy, 19G Soft Closed Tip Catheter (sterile item code 0421157-48)

Combined Spinal Epidural (CSE) Trays

Components

EpiSpin Locking System consists of Tuohy needle with Back-Eye and metal stylet, wings, and locking mechanism along with Sprotte needle with magnifying hub, Clamping Adapter, Epidural Filter, Epidural Label, 25Gx1.5" Hypodermic Needle, 18Gx1.5" Hypodermic Needle, 4" Filterstraw, Glass LOR 5ml Luer Lock Syringe OR Plastic LOR 10ml Luer slip, 3ml Luer Lock Syringe, 5ml Luer Slip Syringe, 20ml Luer Lock Syringe, 3/4 fl. oz. Povidone-Iodine, 3 Prep Sponges, 4"x4.5" Film Dressing, 24"x32" Clear Drape with 4" Circular Opening, 3 Gauze Pads, Paper Towel, CSR Wrap, Needle Foam Block, Drug ID Label, 5ml 1% Lidocaine HCl, 5ml Lidocaine+Epinephrine (Test Dose), 10ml Sodium Chloride, 2ml 0.75% Bupivacaine HCl with 8.25% Dextrose OR Marcaine, where indicated

Luer Item No.	NRFit Item No.	Tray Description
TCL101	TCN101	EpiSpin Lock with 17G x 90mm (3.5") Tuohy needle with 25G x 152mm (6") SPROTTE needle and EpiLong Soft 19G x 90cm Closed Tip catheter (sterile item code 0011351-20)
TCLO201	TCNO201	EpiSpin Lock with 17G x 90mm (3.5") Tuohy needle with 25G x 152mm (6") SPROTTE needle and EpiLong Soft 19G x 90cm Closed Tip catheter (sterile item code 0011351-20)



PAJUNK[®] Choice

We're proud to offer the highest quality Pajunk components, manufactured in Germany and Procedure Trays assembled in the United States.



Also available in **NRFIT**[®]



Custom Trays

Available for **Regional Anesthesia**
and Pain Management Procedures

Pajunk provides a full range of customizable trays where you choose which components you want and how many. Start your Custom Tray today by contacting a Pajunk sales representative.

- PNB
- CPNB
- Support
- Neuraxial
- Chronic Pain



EDUCATIONAL
WEBINARS

ENSURING A HARMONIC PROCEDURE

Regional Anesthesia Procedure Trays



REQUEST SAMPLES NOW!
pajunktrays.com · info@pajunk-usa.com

PAJUNK®

The logo for PAJUNK, featuring the word "PAJUNK" in a bold, teal, sans-serif font. A registered trademark symbol (®) is located to the upper right of the letter "K".

PAJUNK®

PAJUNK® Medical System L.P.
United States of America

4575 Marconi Drive
Alpharetta, GA 30005, USA
Tel.: +1 770.493.6832
Fax: +1 678.514.3388
info@pajunk-usa.com
pajunkusa.com

pajunkusa.com