


**PAJUNK®**

*Safersonic  
Ultrasound Probe Covers with Adhesive*



**Nerve Blocks**

# Safersonic Ultrasound Probe Covers

## Sterile Gel-Free Covers with Adhesive

Thanks to an integrated adhesive strip that provides optimal permeability of the ultrasound waves, the Safersonic probe covers can be attached to the ultrasound probe without the use of any type of contact media (e.g. sterile gel). This reduces the setup and clean up time, and reduces the chance of potential transducer damage due to repeated cleaning. Additionally, the integrated adhesive strip allows air bubble and wrinkle free application of the cover to the probe which results in an improved ultrasound image quality.



### Easy, fast and clean setup

No need for ultrasound gel inside the probe cover due to the integrated adhesive strip



### Easy removal without sticky residues

The probe remains clean throughout the whole procedure resulting in no additional time required for clean up

## Studies

Marhofer P. et al. Sterile working in ultrasonography: the use of dedicated ultrasound covers and sterile ultrasound gel, Expert Rev Medical Devices 2015; 12(6):667-73

Marhofer P. et al. The use of ultrasound-guided regional anaesthetic techniques in neonates and young infants, Acta Anaesthesiol Scand 2014; 58(9):1049-60

Marhofer P. et al. Hygieneaspekte in der ultraschallgestützten Regionalanästhesie, Der Anaesthesist 2016; 65:492-498

Maecten T. Materials Used for Peripheral Nerve Blocks and in Ultrasound Guided Regional Anesthesia, RAPM 2011; 36(7): E6-E8

Product	Item No.	PU
Safersonic Ultrasound Probe Cover with Sterile Gel Pack, 18 x 30cm	S11101	50
Safersonic Ultrasound Probe Cover, 18 x 30cm	S11103	50
Safersonic Ultrasound Probe Cover with Sterile Gel Pack, 18 x 120cm	S11102	50
Safersonic Ultrasound Probe Cover, 18 x 120cm	S11104	50